



Mohawk College
Kickoff

Agenda

Meet the Team	03
SimX Overview + Products	04
Hardware and Software Setup	08
Training Session	11
Timeline	14
Next Steps	15

Meet the Team



**Sarah Friend DNP, ARNP,
FNP-C**

Client Solutions

sarah.friend@simxar.com

Responsible for facilitating
future purchases and
communicating new features



Fedak Arashi

Sales Operations

fedak.arashi@simxar.com

Responsible for project
facilitation and escalation



Leigh Polson

Customer Success



leigh.polson@simxar.com

Responsible for training and
ongoing support


SimX Overview + Products



SimX uses VR to save lives by making simulation training more affordable, easily accessible, and more realistic.



Improves Realism



- Any patient
- Any tools
- Any environment

Increases Flexibility



- Sets up in <5 minutes
- Can be run anywhere
- Multiplayer in person or across the world

Reduces Costs



- 1/10th the cost of mannequin sim
- Automated cases reduce staff workload


Summary of Products

4 Licenses + Standard Curriculum



Undifferentiated Chest Pain (ACLS)
Pediatric Resuscitation (PALS)
Neonatal Resuscitation (NRP)
Motor Vehicle Accident (ATLS)
Shortness of Breath (COVID)
Septic Shock
Preeclampsia
Sandbox

5 Marketplace Cases



Hardware and Software Setup

Technical Setup

134635

Set up your SimX Account on the Oculus Quest


Below is a step-by-step guide on how to set up and run SimX.

In case you don't have one set up, here are instructions for creating a Meta Account.

First, go to <https://auth.meta.com> and click "Set up account with email" and follow the prompts.



After creating a Meta account, please download and log into the Oculus mobile app to finish configuring the headset.




Turn on the headset and follow the prompts. When configuring the headset, you will be prompted to create a Meta Horizon Profile. Feel free to skip creating an avatar, payment, etc.

Step 2: Install the SimX App on the headset(s):

Visit [HERE](#) to find the SimX app in the app store on your computer.

Step 3: Install the Moderator Client on your Windows Computer:

You should have received a message from SimX to your institutional email inviting you to create a SimX account. Follow the link to complete your registration.



Network Setup


To ensure the training session will run smoothly, please share the following documents with your IT department.



End User Network Connection Test

After it has been confirmed that your IT department has whitelisted the domains provided in the SimX System Requirements PDF, each SimX Moderator Client has the ability to test the network connection. To do so, follow the instructions below:

1. When in the Moderator Client, select the green gear icon located at the top left corner of the screen.



SimX End User Network
Connection Test



Protocols/Websites/Network Ports/IPs

SimX VR headset and computer applications communicate via HTTPS websockets with the SimX cloud servers to download assets and synchronize scenario status.

For customers operating on secured networks, it may be necessary to explicitly set the network security systems to allow the catch-all root domain *simx-infra.com*.

For some systems, it may be necessary to explicitly allow all of the individual SimX server URLs. A comprehensive list of these URLs is kept updated here: <https://static-assets.simx-infra.com/simx-allow-list-full.txt>

Bandwidth

We don't generally provide detailed requirements because SimX is very bandwidth insensitive as content is loaded onto devices and rendered on-device with minimal communication with our central server. SimX only requires external connectivity over standard web ports (80 and 443).

The one instance where low bandwidth can cause a sub-optimal experience is during initial loading of content onto headsets, as some are a few hundred mb per scenario.

However, this only needs to occur once with every version of the scenario and then the content is permanently cached on the headset and so we recommend that content for any given session be preloaded onto the headset.

SimX System
Requirements

Training Session

Where:



**What You
Need:**



Duration:


2.5 Hrs

Headset Setup




During this period, we will work with the POCs to set up their Play Space and sync Multiplayer on the headset.

Moderator Setup



During this period, we will work with the POCs to launch the SimX Moderator Application on the laptop.

Timeline

	4/11	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10
Kickoff & Setup	Kick Off	Create Accounts / Setup Headsets				Ready for SimX Cases					
IT / Network Setup		Whitelist Domains									
Training				Virtual Training Session							

Next Steps:

Mohawk College

- Setup headsets and download SimX
- Schedule virtual training

SimX

- Send presentation
- Send setup instructions
- Send virtual training Calendly Link



SIM^x

Thank you!